



cmd_key™

Programmed for QLab 4

Congratulations on your **cmd_key**!

The **cmd_key** for QLab keyboard is pre-programmed from the factory to work with QLab 4 from Figure 53. It is programmed using the default keyboard shortcuts provided with QLab 4.

You can assign Hot Keys 1-5 and A-F for cue triggers or scripting of your choice. Please refer the the QLab manual on how to assign them.

Hot Key			Hot Key		
1	=	!	A	=	^
2	=	@	B	=	&
3	=	#	C	=	*
4	=	\$	D	=	(
5	=	%	E	=)
			F	=	_

When you plug the keyboard into the Mac for the first time, MacOS might ask you to identify the keyboard by pressing a specific key on the keyboard. Just close that window and your Mac will treat it as a QWERTY keyboard, which is how it's programmed.

cmd_key is programmed for QWERTY US Layout only.

The hub underneath the **cmd_key** can be used for a Mouse. It is not designed for a USB Thumb Drive. It is underpowered for that type of device.

If you need tech support, please reach out to us.

www.cmd-key.com

www.highoutput.com

cmd_key@highoutput.com

Figure 53® and QLab® are registered trademarks of Figure 53, LLC. cmd_key and High Output, Inc. are not affiliated with Figure 53, LLC and this keyboard has not been reviewed nor is it approved by Figure 53, LLC